




# SOP 005 - Preparation of Tissue Homogenates

## Table of Contents

1. Principle.....	3
2. Specimens.....	3
3. Reagents and Materials.....	3
4. Standards, Controls, and Solutions.....	3
5. Equipment and Special Supplies .....	3
6. Procedure.....	3
7. References .....	4

Uncontrolled Copy

## SOP 005 - Preparation of Tissue Homogenates

SOP Name: <b>Preparation of Tissue Homogenates</b>		SOP #: <b>005</b>
North Carolina Office of the Chief Medical Examiner Toxicology Laboratory	<b>Revision:</b>	<b>Revision Date/Initials:</b>
<b>Approving Authority Name</b>	<b>Approving Authority Signature</b>	<b>Approval Date</b>
Ruth E. Winecker, Ph.D.		04/07/2015
Ruth E. Winecker, Ph.D.		06/06/2016
Ruth E. Winecker, Ph.D.		08/21/2017

## SOP 005 - Preparation of Tissue Homogenates

### 1. Principle

- 1.1. procedure describes the method for preparing 1:4 tissue homogenates for determining drug concentrations in the tissues. Tissues obtained from the pathologist or medical examiner are diluted with water and completely homogenized.

### 2. Specimens

- 2.1. This procedure is applicable to liver, brain, spleen, kidney, muscle, or other tissue specimens.

### 3. Reagents and Materials

- 3.1. Deionized water
- 3.2. Tissue
- 3.3. Disposable scalpel
- 3.4. 50mL Falcon Tubes
- 3.5. 10% Bleach Solution

### 4. Standards, Controls, and Solutions

- 4.1. N/A

### 5. Equipment and Special Supplies

- 5.1. Top-loading balance
- 5.2. Silent Crusher M homogenizer

### 6. Procedure

- 6.1. Note: When preparing an hd2 (second homogenate prepared from parent liver specimen) clean homogenizer generator thoroughly (6.10) prior to next step.
- 6.2. In a 50mL Falcon Tube, weigh approximately 4g of tissue on the top-loading balance using a scalpel to cut it into small enough pieces to fit in the generator (Approx. 1/8 inch).
- 6.3. Without taring the balance, multiply the exact amount of tissue weighed by four.

## SOP 005 - Preparation of Tissue Homogenates

- 6.3.1. Add deionized water to the Falcon tube containing the tissue until the final mass reaches the calculated mass.
  - 6.3.2. **E.G. 4.0g tissue (4g x 4 = 16g) → add water to a final mass of 16g.**
  - 6.4. Place the Falcon tube under the generator of the Silent Crusher M homogenizer.
  - 6.5. Lower the generator into the Falcon tube so that it rests approximately ¼ inch off the bottom and tighten the adjustment knob to hold it in place.
  - 6.6. Turn the Silent Crusher M on, press the start button, and increase the speed to 10000 RPM.
  - 6.7. Blend the tissue until completely homogenized. (1-2 min)
  - 6.8. Pour the contents of the bottle into a 20 mL scintillation vial labeled with T#, S#, contents and dilution factor (e.g. 1:4 hd1 liver), date prepared, and initials of preparer.
  - 6.9. Store in the refrigerator at about 4°C.
  - 6.10. Clean the homogenizer by disassembling the generator and cleaning all parts thoroughly with bleach solution and clean running DI water.
7. **References**
- 7.1. Levine, Barry. "Postmortem Forensic Toxicology." *Principles of Forensic Toxicology*. 2nd ed. Washington, DC: AACC, 2006. 3-13. Print.
  - 7.2. Heildolph Instruments, comp. *SilentCrusher M*. 14 Oct. 2011. Instruction manual. Room 2606 OCME Toxicology, Raleigh, NC.