

Technical Procedure for the CrimeScope and Mini-CrimeScope

1.0 Purpose - This procedure describes how to examine evidence with the CrimeScope and Mini-CrimeScope alternate light sources.

2.0 Scope - This procedure applies to all evidence that is examined with the CrimeScope and Mini-CrimeScope.

2.1 The CrimeScope and Mini-CrimeScope can assist in searching for latent fingerprints, palmprints, footprints, footwear impressions and other impressions. Fingerprint residue, certain biological fluids, trace evidence and other compounds and substances may fluoresce naturally when exposed to different wavelengths of light. In addition, the CrimeScope and Mini-CrimeScope may be used with fluorescent dyes to develop and enhance latent impressions. The CrimeScope and Mini-CrimeScope may be used in the laboratory and when providing technical field assistance.

3.0 Definitions – N/A

4.0 Equipment, Materials and Reagents

4.1 Equipment and Materials

- CrimeScope, Mini-CrimeScope, and protective goggles
- Camera equipment

4.2 Reagents – N/A

5.0 Procedure – Items of evidence shall be subjected to the light source prior to using processing techniques, if utilizing the light source for inherent luminescence examination. This will serve to detect any inherent latent impressions and to reveal the color and the intensity of any background fluorescence. The background fluorescence will allow the analyst to determine which fluorescent dye will be appropriate to use in the sequence of processing evidence. The CrimeScope and/or Mini-CrimeScope shall also be used in conjunction with fluorescent dyes and chemicals to develop and photograph latent impressions.

5.1 Start Up Procedures

5.1.1 Plug the unit into a standard electrical outlet.

5.1.2 Attach the fiber optic cord.

5.1.3 Turn the power rocker switch ON. This will activate the fan system of the unit (the fan sound will be evident in a few seconds).

5.1.4 Turn ON the Power Control Switch. This will activate the lamp (a clicking sound will indicate the lamp is activated and the unit will initiate the light beam).

Note: The Mini-CrimeScope only has a power control switch.

5.2 Shut Down Procedures

5.2.1 Turn the Power Control Switch to the OFF position.

5.2.2 Allow the unit to cool, with the fan running, for approximately 10-15 minutes.

5.2.3 After the unit has cooled, check the body of the unit to ensure that it is cool.

5.2.4 Turn the fan power switch to the OFF position.

Note: The Mini-CrimeScope only has a Power Control switch and does not need the additional run time for the fan.

5.3 Wavelength and Filter Selection

5.3.1 To select the desired wavelength, turn the wavelength control or selector to the desired setting.

5.3.2 The wavelength in use will be displayed on the front of the unit and/or on the optic handheld attachment.

5.3.3 Various wavelengths and filter selections shall be used to determine which is most appropriate for the item of evidence being examined.

5.4 Examination of Evidence

5.4.1 While wearing goggles, scan the item of evidence with the light source. A strong hand held magnifier may be utilized to enhance visualization of the latent impressions.

5.4.2 If a latent impression is developed, immediately note the location on the item and direct the light source away from the area.

5.4.3 Position the area of interest under the camera and place the appropriate filter over the lens of the camera (see technical procedures for the camera used).

5.4.4 Direct the light source over the area to be photographed and begin photography.

5.4.5 After each photograph is taken, redirect the source away from the impression to avoid destruction of the area or surface. Prolonged exposure to the beam may cause the latent impression to photo-degenerate over a short period of time until eventually the impression disappears.

5.4.6 After the impression is photographed, further techniques may be applied to enhance or develop additional impressions, if applicable. (See Technical Procedures for Friction Ridge Analysis and Comparison)

5.5 Standards and Controls – N/A

5.6 Calibration – N/A

5.7 Sampling – N/A

5.8 Calculations – N/A

5.9 Uncertainty of Measurement – N/A

6.0 Limitations - N/A

7.0 Safety – Eye protection shall be worn when utilizing a laser and/or alternate light source. This includes any individual who may be in the same room or area. This is particularly important when examining reflective surfaces as the light source may reflect and result in eye damage.

7.1 Never look directly into any light source as this will cause eye damage.

7.2 Do not expose the light source to the skin; it will not immediately cause harm, but prolonged exposure may have long term effects.

8.0 References

Almog, J., A. Hirshfeld and J.T. Klug. "Reagents for the Chemical Development of Latent fingerprints: Synthesis and Properties of Some Ninhydrin Analogues." *Journal of Forensic Sciences*. Vol. 27, 4: 912-917 (1982).

Ayala, Garcia J.A. "Lasers in Crime Investigation: Detecting Latent Fingerprints at the Scene of a Crime." *International Criminal Police Review*. (April 1984): 90-93.

Burt, J.A. and E.R. Menzel. "Laser Detection of Latent Fingerprints: Difficult Surfaces." *Journal of Forensic Sciences*. Vol. 13, 2: 364-370 (1985).

Dalrymple, B.E., J.M. Duff and E.R. Menzel. "Inherent Fingerprint Luminescence – Detection by Laser." *Journal of Forensic Sciences*. Vol. 22, 1: 106-111 (1977).

Everse, K.E. and E.R. Menzel. "Blood Print Detection by Fluorescence." *Fluorescence Detection*. Center for Forensic Studies, Texas Tech University, Lubbock. TXSPIE. Vol. 743: 184-202 (1987).

Fisher, J.F. "The Enhancement of Blood Prints by Chemical Methods and Laser-Induced Fluorescence." *Identification News*. (1984): 2; 14-15.

Forensic Laser Technology. Omnichrome. New Jersey, 1-6.

German, E.R. "You Are Missing Ninhydrin Developed Prints." *Identification News*. (1981): 3-5.

Goss, E., L. Sin-David and J. Almog. "Transmitted Infrared Luminescence in Document Examination." *Journal of Forensic Sciences*. Vol. 25, 2: 382-385 (1980).

Hammond, J. "Cyanoacrylate Ester Fuming For the Development of Latent Prints." *Loctite Corporation*: 1-24.

Hazen, R.J. "Significant Advances in the Science of Fingerprints." *Forensic Science*. 2nd Ed. Geoffrey Davies, Editor ACS, Washington, DC. (1986): 299-312.

Herod, D.W. and E.R. Menzel. "Laser Detection of Latent Fingerprints: Ninhydrin Followed by Zinc Chloride." *Journal of Forensic Science*. Vol. 27, 3: 513-518 (1982).

Inlow, V.K. "The Use of Flashlight and Filters for Fluorescent Examinations." *Administrative Advanced Latent Fingerprint Training, FBI, Quantico, VA*. (November 1993): 1-5.

“III: Visualization of Latent Fingerprints by Fluorescent Reagents in Vapor Phase.” *Journal of Forensic Sciences*. Vol. 25, 2: 408-410 (1980).

Kobus, Hilton J., Ronald N. Warrenner and Milutin Stoilovic. “Laser Detection of Latent Fingerprints: Treatment with Glue Containing Cyanoacrylate Ester.” *Journal of Forensic Sciences*. Vol. 28, 2: 307-317 (1983).

“Laser Enhancement of Latent Fingerprints.” Fingermatrix, Inc., Pender Drive, Suite 310 Fairfax, Virginia, 22030. (March 17, 1988): 1-16.

“Luminescent Enhancement Procedures for the Detection of Latent Fingerprints.” *Proceedings of the International Forensic Symposium on Latent Prints*. (July 1987): 45-49.

Menzel, E.R. “Detection of Latent Fingerprints by LASER-Excited Luminescence.” *Analytical Chemistry*. Vol. 61, 8: 557 – 561 (1989).

Menzel, E.R. “Laser Fingerprint Detection and Development.” *Proceedings of the International Forensic Symposium on Latent Prints*. (July 7-10, 1987): 25-38.

Menzel, E.R. “Pretreatment of Latent Prints for Laser Development.” *Forensic Science Review*. Vol. 1, 1: 43-66 (June 1989).

Menzel, E.R. “A Guide to LASER Latent Fingerprint Development Procedures.” *Identification News*. (September 1983): 9 – 13.

Menzel, E.R. “LASER Detection of Latent Fingerprints – Treatment with Phosphorescers.” *Journal of Forensic Sciences*. Vol. 24, 3: 582 – 585 (1979).

Menzel, E.R. Chemical Reagents for the Development of Latent Fingerprints.

Menzel, E.R., and J.M. Duff. “Laser Detection of Latent Fingerprints – Treatment with Fluorescers.” *Journal of Forensic Sciences*. Vol. 24, 1: 96-100 (1979).

Menzel, E.R. “Ion-Laser Detection of Fingerprints Grows More Powerful.” *Laser Focus World*. (November 1989): 89-94.

Menzel, E.R. “Laser Detection of Latent Fingerprints on Skin.” *Journal of Forensic Sciences*. Vol. 27, 4: 918-922 (1982).

Menzel, E.R., et al. “Room Light and Laser Development of Latent Fingerprints with Enzymes.” *Journal of Forensic Sciences*. Vol. 29, 1: 99-109 (January 1984).

Menzel, E.R. and J.A. Burt. “LASER Detection of Latent Fingerprints: Difficult Surfaces.” *Journal of Forensic Sciences*. Vol. 13, 2: 364 – 370 (1985).

Misner, A.H. “Ultraviolet Light Sources and Their Uses.” *Journal of Forensic Identification*. Vol. 41, 3: 171-175 (June 1991).

Salares, V. R., C.R. Eves and P.R. Carey. “On the Detection of Fingerprints by Laser Excited Luminescence.” *Forensic Science International*. Vol. 14: 229-237 (1979).

Sasson, Y. and J. Almog. "Chemical Reagents for the Development of Latent Fingerprints, I: Scope and Limitations of the Reagent 4-Dimethylamino-Cinnamaldehyde." *Journal of Forensic Sciences*. (1978): 852-855.

Thompson, J. *Laser Technology as Utilized in the Development of Latent Prints*. Ocean County Sheriff's Department, Criminalistics Investigation Unit.

Walton, A. N. "A Technique for the Detection and Enhancement of Latent Prints on Curved Surfaces by the Use of Fluorescent Dyes and Painting with Laser Light (Beam)." *Proceedings of the International Forensic Symposium on Latent Prints*. (July 7-10, 1987): 121-123.

Watkin, J.E. "Alternate Lighting Methods of Detecting Latent Prints." *Proceedings of the International Forensic Symposium on Latent Prints*. (July 7-10, 1987): 39-44.

9.0 Records – N/A

10.0 Attachments – N/A

Revision History		
Effective Date	Version Number	Reason
09/17/2012	1	Original Document
10/31/2013	2	Added issuing authority to header
04/07/2017	3	Header Update – Removed Digital reference. 5.4.6 – updated procedure reference.
02/01/2019	4	Added "Mini-CrimeScope" in addition to anywhere "CrimeScope" was mentioned throughout document 1.0: spelling correction 5.1.4: added note; lowercase "power control switch" 5.2.1: changed "lamp power switch" to "power control switch" 5.2.2: changed numbers to numerals instead of spelling out 5.2.4: added note