




SOP 071 Routine LC-MS Maintenance

Table of Contents

1. Purpose	3
2. Reagents and Materials.....	3
3. Procedure for Ion Transfer Tube and Sweep Cone Cleaning.....	3
4. Procedure for Trimming ESI Fused-Silica Sample Tube.....	4
5. References	4

Uncontrolled Copy

SOP 071 Routine LC-MS Maintenance

SOP Name: Routine LC MS Maintenance		SOP #: 071
North Carolina Office of the Chief Medical Examiner Toxicology Laboratory	Revision:	Revision Date/Initials:
	3.5 – Deleted 4.1, 5.1, 5.3 – Updated Hyperlinks 5.2 – Added mention of hard copy maintenance procedures.	JOB 09/12/2017 MSF 09/12/2017
Approving Authority Name	Approving Authority Signature	Approval Date
Ruth E. Winecker, Ph.D.		04/15/2015
Ruth E. Winecker, Ph.D.		06/10/2016
Ruth E. Winecker, Ph.D.		09/12/2017

SOP 071 Routine LC-MS Maintenance

1. Purpose

- 1.1. This SOP gives the procedures and references for routine maintenance of Thermo Scientific LC/MS analytical instruments.

2. Reagents and Materials

- 2.1. Nitric Acid
- 2.2. Distilled water (HPLC grade or 18.2 M Ω ·cm preferred)
- 2.3. HPLC grade methanol
- 2.4. HPLC grade acetone (optional)
- 2.5. 50mL graduated cylinder
- 2.6. Paper towels and/or Kimwipes
- 2.7. Hypodermic needle (P/N 00106-20000 RevB, 28 gauge RW)
- 2.8. Ion transfer tube removal tool (P/N 70111-20258)

3. Procedure for Ion Transfer Tube and Sweep Cone Cleaning

- 3.1. Turn instrument to standby
 - 3.1.1. Note: It is not necessary to vent the system or reduce the capillary temperature.
- 3.2. Remove ion source housing from front of the MS detector (see respective instrument hardware manual for more details).
- 3.3. Remove the ion sweep cone by grabbing the outer ridges and pulling it straight off. Set aside to cool.
 - 3.3.1. Caution: the sweep cone is very hot; use paper towels to help grab the sweep cone to avoid getting burned
 - 3.3.2. Clean by rinsing (can gently rub using paper towels or Kimwipes) with distilled water followed by methanol.
- 3.4. Using the special removal tool, remove ion transfer tube by turning it counterclockwise until you can pull it free. Set aside to let cool for 5 minutes.
 - 3.4.1. Be careful not to lose the O-ring that fits between the ion transfer tube and the spray cone.
- 3.5. Sonicate the ion transfer tube for 30-45 minutes in a solution of nitric acid (15-20% v/v) (a 50mL graduated cylinder works well for this step).
- 3.6. Rinse with distilled water.
- 3.7. Sonicate for 5-10 minutes in methanol.

SOP 071 Routine LC-MS Maintenance

- 3.7.1. Use HPLC grade methanol (or better) that has not come into contact with plastics such as the repeater pipette tips.
 - 3.8. (Optional) Sonicate for 5 minutes in acetone.
 - 3.8.1. Use HPLC grade (or better) that has not come into contact with plastics.
 - 3.9. Replace ion transfer tube (ensure that the O-ring is present), carefully rotating as it is inserted.
 - 3.10. Using the special removal tool, rotate clockwise () until it is finger tight.
 - 3.11. Reinstall the ion sweep cone.
 - 3.12. Reinstall the ion source housing.
4. **Procedure for Trimming ESI Fused-Silica Sample Tube**
 - 4.1. See S:\toxicology\QAQC\SOP\Media\IonMax_API_IonSouce_Hardware.pdf (page 12-13).
5. **References**
 - 5.1. [LXQ Hardware Manual](#)
 - 5.2. Hard copies of maintenance procedures can be found in binders located next to each instrument.
 - 5.3. Ion Max and Ion Max-S API Source Hardware Manual
 - 5.3.1. S:\toxicology\QAQC\SOP\Media\IonMax_API_IonSouce_Hardware.pdf