




SOP 100 - Salicylate Colorimetric Screen

Table of Contents

1. Principle of Assay.....	3
2. Specimens.....	3
3. Reagents and Materials.....	3
4. Standards, Controls, and Solutions.....	3
5. Equipment and Special Supplies	4
6. Procedure.....	5
7. Calculations	5
8. Quality Control.....	5
9. Load Assignment Packet Preparation.....	5
10. References.....	6

Uncontrolled Copy

SOP 100 - Salicylate Colorimetric Screen

SOP Name: <h3 style="text-align: center; margin: 0;">Salicylate Colorimetric Screen</h3>	SOP #: <h3 style="text-align: center; margin: 0;">SOP 100</h3>	
North Carolina Office of the Chief Medical Examiner Toxicology Laboratory	Revision:	Revision Date/Initials:
Approving Authority Name	Approving Authority Signature	Approval Date
Ruth E. Winecker, Ph.D.		04/08/2015
Ruth E. Winecker, Ph.D.		06/10/2016
Ruth E. Winecker, Ph.D.		08/29/2017

SOP 100 - Salicylate Colorimetric Screen

1. Principle of Assay

- 1.1. This method is designed to detect the presence of salicylate in blood, urine, or tissue specimens by way of a colorimetric test. The specimens are prepared with Trinders reagent and a purple color indicates a positive result.
- 1.2. The salicylates are a group of compounds used for the purposes of analgesia, fever reduction, and decreasing inflammation. In the blood acetyl salicylic acid (ASA) is rapidly hydrolyzed to salicylic acid (SA). Salicylate has a half-life that is dose-dependent and variable, ranging from 2-3 hours to 20-30 hours. In overdose, the half-life is greatly increased because the metabolic acidosis causes both a decrease in both blood and urine pH which reduces salicylate elimination.

2. Specimens

- 2.1. Blood, vitreous, serum, Plasma, urine, bile, gastric contents, or tissue homogenate - volume varies depending upon the specimen availability, typically 1 mL of blood, urine, vitreous, bile, or gastric contents, or 1 g of tissue homogenate (1:4 dilution).

Specimen	Aliquot Volume
Blood	1.0 mL
Vitreous Humor	1.0 mL
Serum	1.0 mL
Plasma	1.0 mL
Urine/Bile/Gastric	1.0 mL
Tissue homogenate	1.0 g

3. Reagents and Materials

- 3.1. Mercuric chloride
- 3.2. Deionized water
- 3.3. Ferric nitrate
- 3.4. Salicylate standard (1000 mg/dL)
- 3.5. Blood bank blood (drug-free).
- 3.6. Hydrochloric Acid

4. Standards, Controls, and Solutions

SOP 100 - Salicylate Colorimetric Screen

4.1. Salicylate Standard (1000mg/dL)

- 4.1.1. Weigh 117 mg of sodium salicylate using an analytical balance.
- 4.1.2. Place in a 10 mL volumetric flask and dilute to the mark with deionized water. Mix Well

4.2. 1N Hydrochloric Acid

- 4.2.1. In a 250 mL volumetric flask, place about 100 mL of deionized water.
- 4.2.2. Slowly add 21 mL of concentrated hydrochloric acid and gently vortex.
- 4.2.3. Dilute to the mark with deionized water.
- 4.2.4. Invert several times to mix.

4.3. Trinders reagent

- 4.3.1. In a 2000 mL Erlenmeyer flask with a stir bar, place 40 g of mercuric chloride, and 850 mL of deionized water.
- 4.3.2. Heat and stir until warm to the touch.
- 4.3.3. Cool and add 120 mL of 1 N HCl and 40 g ferric nitrate.
- 4.3.4. Pour into a 1000 mL volumetric flask and dilute to the mark with deionized water. Mix well.

5. Equipment and Special Supplies

- 5.1. 5 mL Disposable Centrifuge Tubes
- 5.2. Vortex mixer
- 5.3. Centrifuge
- 5.4. Micro Pipette, air or positive displacement 10-100 μ L, calibrated
- 5.5. Micro Pipette, air or displacement 100-1000 μ L, calibrated
- 5.6. Macro Pipette, air displacement 1-5mL

SOP 100 - Salicylate Colorimetric Screen

6. Procedure

- 6.1. Label four 5mL Disposable Centrifuge Tubes "20", "50", "blank", and "QC Std". Do the same for all assigned specimens. Pipette 1mL of blank blood into the four control tubes.
- 6.2. Add 20uL of 1000mg/dL salicylate standard to tube labeled "20".
Add 50uL of 1000mg/dL salicylate standard to tube labeled "50".
Add 20uL of 1000mg/dL salicylate QC standard to tube labeled "QC Std".
- 6.3. Pipette 1mL of specimen (unless otherwise directed upon assignment) into appropriately labeled 5mL Disposable Centrifuge Tubes.
- 6.4. Add 5mL Trinders Reagent to all standards and specimens, cap, vigorously shake for 30 seconds, and centrifuge at 2000 x g for 10 minutes.
- 6.5. Remove from centrifuge and observe.
- 6.6. A positive result is indicated by a purple color change. A negative result is indicated by no change in Trinders Reagent (yellow).
- 6.7. Have another analyst or toxicologist observe the results and initial the result sheet

7. Calculations

- 7.1. N/A

8. Quality Control

- 8.1. For an analysis to be acceptable the following criteria must be met:
 - 8.1.1. A positive result (purple color change) must be observed in the 20 and 50mg/dL standards as well as the 20mg/dL QC.
 - 8.1.2. A negative result (no color change to the yellow Trinders reagent) must be observed in the negative control.
 - 8.1.3. If the criteria above are not met, the specimens will be re-analyzed.

9. Load Assignment Packet Preparation

- 9.1. After completion of the analysis, the analyst will organize the data in the following order:

SOP 100 - Salicylate Colorimetric Screen

- 9.1.1. The Load Checklist - should be initialed and dated to acknowledge completion of load.
- 9.1.2. Load assignment sheets, followed by any additional notes to file pertaining to load.
- 9.1.3. Load specimen sheet.
- 9.1.4. Chain of Custody.
- 9.1.5. Aliquot Chain of Custody.
- 9.1.6. Standard and Control Worksheet.
- 9.1.7. Results worksheet – a copy for each case in the load and one for the Tox Folder.

10. References

- 10.1. Levine, Barry. *Principles of Forensic Toxicology*. Washington, D.C.: American Association for Clinical Chemistry, 1999. Print.
- 10.2. Clarke's Analysis of Drugs and Poisons. 3rd Edition. Moffat, Osselton and Widdop eds. Pharmaceutical Press, London. 2004. P. 296.