



**SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for
Quantification by Liquid Chromatography/Electrospray Mass
Spectrometry/ Mass Spectrometry (LC/MS/MS)**

Table of Contents

1. Principle of Assay.....	3
2. Specimens.....	3
3. Reagents and Materials.....	4
4. Standards, Controls, and Solutions.....	4
5. Equipment and Special Supplies	5
6. Instrumentation and Parameters	5
7. Target Ions.....	5
8. Procedure.....	5
9. Calculations	6
10. Quality Control	8
11. Validation of Method.....	9
12. Reporting.....	10
13. Preparation of Load.....	11
14. References.....	12

**SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for
Quantification by Liquid Chromatography/Electrospray Mass
Spectrometry/ Mass Spectrometry (LC/MS/MS)**

SOP Name: Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)		SOP #: 416
North Carolina Office of the Chief Medical Examiner Toxicology Laboratory	Revision:	Revision Date/Initials:
	Initial document preparation 11. – Added Validation table 12.4.1 – Updated reporting units 12.4.1 – Updated reporting units (ng/mL) 4 – Added standard prep instructions 6 – Updated instrument parameters 9.1.1 – Updated instrument method name	MSF – 05/10/2016 MSF – 06/20/2016 MSF – 11/21/2016 MSF – 05/25/2017 MSF – 06/05/2017 MSF – 08/25/2017
Approving Authority Name	Approving Authority Signature	Approval Date
Ruth E. Winecker, Ph.D.		08/25/2016
Ruth E. Winecker, Ph.D.		12/07/2017

SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)

1. Principle of Assay

- 1.1. This method is designed to confirm and quantitate furanylfentanyl in biological specimens by Liquid Chromatography Tandem Electrospray Mass Spectrometry (LC-MS/MS). Furanylfentanyl is extracted from biological matrices by protein precipitation with acetone and identified by retention time and ion ratio of product ions.
- 1.2. Furanylfentanyl is another drug in the ever-growing class of “novel psychoactive substances” (NPS, note: this abbreviation is often used, synonymously, with “new psychoactive substances”). The vast majority of NPS drugs are derivatives of known legal or illicit drugs, and are meant to circumvent existing drug laws. They are commonly sold in retail stores or on-line, often with the caveat of “not meant for human consumption.” The design and synthesis of such drugs are pulled from patents and scientific publications, with very little known about activity or toxicity in humans.

First identified in the United States in the latter half of 2015, furanylfentanyl shares the same structural backbone of fentanyl. There are no known pharmacokinetic, pharmacodynamic, or toxicity studies in higher mammals or humans. Like fentanyl, however, furanylfentanyl is expected to act as an opioid and have central analgesic effects, as well as adverse effects of respiratory depression.

Because it has not been studied in humans, there are no known “therapeutic doses”, and due to the way it is produced and sold, users may not know the concentration or purity of furanylfentanyl, nor even know they are using it at all, if it is being sold as an adulterant to heroin or claimed to be another drug. Outside of this laboratory, two deaths in Sweden have been associated with furanylfentanyl intoxication, with serum concentrations of 4.4 and 148 ng/mL (1).

Like fentanyl, furanylfentanyl is expected to exhibit postmortem redistribution (2, 3). As with any drug, particularly one commonly likely to be combined with other opioid analgesics, and CNS depressants such as benzodiazepines, interpretation of postmortem furanylfentanyl concentrations must rely not only on the drug concentration, but incorporate patient history, autopsy and scene findings.

2. Specimens

- 2.1. This procedure is applicable to urine, blood, serum, *bile, *gastric contents, and properly prepared **tissue specimens (typically 1:4 homogenates). A 0.1 mL (g) specimen amount (in duplicate) is generally employed unless a dilution is required so that the calibration curve encompasses the expected range of unknown specimens.

SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)

Note: Liver is not a viable specimen for analysis due to extensive metabolism/biotransformation of furanylfentanyl *in vivo* and possible *in vitro* stability problems.

- 2.1.1. *For non-typical matrices, an additional 0.1mL aliquot shall be taken (volume permitting), spiked with appropriate QC, and analyzed to help to identify any matrix effects. (See Non-Matched Matrix Protocol section of the QA/QC manual).

3. Reagents and Materials

- 3.1. DI water, HPLC grade
- 3.2. Methanol, HPLC grade
- 3.3. Acetone, HPLC grade
- 3.4. Deuterated Internal Standard Mix
- 3.5. Calibration Standard
- 3.6. QC Standard
- 3.7. Drug Free Blood, Urine, Liver Homogenate
- 3.8. Water with 0.1% formic acid
- 3.9. Methanol with 0.1% formic acid

4. Standards, Controls, and Solutions

4.1. Fentanyl-d5 Stock Solution A (10µg/mL)

- 4.1.1. Into a 10mL volumetric flask, add the contents of 1 ampule (~1mL) of Fentanyl-d5 (Cerilliant – 100µg/mL).
- 4.1.2. Fill to the line with methanol, insert stopper and invert three times to mix. Transfer to properly labeled 16x100mm screw topped test tubes and cap. Store in laboratory refrigerator (R1-2601). See [SOP-010](#).

4.2. Fentanyl-d5 Stock Solution B (1000ng/mL)

- 4.2.1. Into a 10mL volumetric flask, add 1 ml of Fentanyl-d5 Stock Solution A (10µg/mL) with a micropipette.
- 4.2.2. Fill to the line with methanol, insert stopper and invert three times to mix. Transfer to properly labeled 16x100mm screw topped test tubes and cap. Store in laboratory refrigerator (R1-2601). See [SOP-010](#).

4.3. Fentanyl-d5 Internal Standard (10ng/mL)

- 4.3.1. Into a 10mL volumetric flask, add 0.1 ml of Fentanyl-d5 Stock Solution B (1000ng/mL) with a micropipette.

SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)

- 4.3.2. Fill to the line with methanol, insert stopper and invert three times to mix. Transfer to properly labeled 16x100mm screw topped test tubes and cap. Store in laboratory refrigerator (R1-2601). See [SOP-010](#).
- 4.4. **Furanylfentanyl Calibrators and Positive Controls** – these standards are to be prepared by the QA/QC Chemist or appointee. Inform the QA/QC Chemist if calibration/control standards need to be made.
- 4.5. **Water with 0.1% formic acid**
 - 4.5.1. To a 4L bottle of HPLC grade water, add 4 mL of formic acid
 - 4.5.2. Label bottle as “LC/MS” and “with 0.1% formic acid”.
- 4.6. **Methanol with 0.1% formic acid**
 - 4.6.1. To a 4L bottle of HPLC grade methanol, add 4 mL of formic acid
 - 4.6.2. Label bottle as “LC/MS” and “with 0.1% formic acid”.
5. **Equipment and Special Supplies**
 - 5.1. Test Tubes, 16 x 100 mm
 - 5.2. LC autosampler vials, 12 x 32 mm
 - 5.3. Polyspring inserts, 5 mm O.D.
 - 5.4. Centrifuge 2000 x g
 - 5.5. Vortex mixer
 - 5.6. Nitrogen evaporator
6. **Instrumentation and Parameters**
 - 6.1. Windows PC with Thermo LCQuan and Xcaliber software
 - 6.1.1. Instrument method (TSQ03): “Furanylfentanyl”
 - 6.1.2. Click [here](#) for instrument parameters.
 - 6.2. Thermo Surveyor LC autosampler, or equivalent
 - 6.3. Thermo Surveyor LC pump, or equivalent
 - 6.4. Thermo TSQ triple quadrupole mass spectrometer
7. **Target Ions** (± 1 nominal mass)
 - 7.1. Furanylfentanyl **(375 188 105)**
 - 7.2. Fentanyl-d5 **(342 188 105)**
 - 7.2.1. Note: The precursor ion of each analyte is listed first and bolded, the first product ion- used for quantification-is second, followed by the second product ion-used for qualification/confirmation.
8. **Procedure**

SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)

- 8.1. Prepare a colored tape label for each standard, blank, control, and specimen to be placed on 16x100 mm test tubes.
- 8.2. Add the appropriate quantity (according to the [Standard & Control Worksheet](#)) of Deuterated Internal Standard Mix to all the tubes.
- 8.3. Add the appropriate quantity (according to the [Standard & Control Worksheet](#)) of calibration standard and QC to the tubes labeled as standards and control, respectively, labeling test tubes as you go. Only internal standard should be present in the test tube labeled “Blank”.
- 8.4. Add 0.1mL of blank blood to all standards, controls, and blank test tubes (0.1 mL blank urine/0.1g blank liver homogenate to urine/liver blank and QC test tubes respectively).
- 8.5. Add the appropriate amount of predetermined unknown specimen labeling test tubes as you go. (See [Specimens](#) section).
- 8.6. Vortex all test tubes for 10 seconds.
- 8.7. Add 3.5mL acetone to each tube and vortex for 20 seconds.
- 8.8. Centrifuge at 2000 x g for 10 minutes.
- 8.9. Decant the top acetone layer into clean and labeled 16x100 test tubes, place in nitrogen evaporator, and evaporate at 55° C to dryness.
- 8.10. Remove dried specimens from nitrogen evaporator and reconstitute with 300µL of methanol. Vortex for 10 seconds and centrifuge at 2000 x g for 5 minutes.
- 8.11. Transfer approximately 100µL of each extract to appropriately labeled autosampler vials fitted with 200 µL polyspring insert and place in the autosampler tray of the Thermo TSQ triple-quadrupole LC/MS/MS.
- 8.12. Build and initiate sequence as directed in [SOP 053](#).

9. Calculations

9.1. Quantification

- 9.1.1. The method for processing the data using the Thermo LCQuan software is “Furanylfentanyl Quant” ([SOP 055](#)). It is used to calculate the internal standard response ratios, raw amounts, concentration, and ion ratios.

SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)

9.1.2. These calculations are computed as follows:

9.1.2.1. Response Ratio:

9.1.2.1.1. Response Ratio = response of the analytes quantifying product ion compared to that of the internal standard's quantifying product ion.

9.1.2.1.2. Response Ratio = QN_a / QN_{istd}

9.1.2.1.2.1. QN_a = response of the quantitative ion of the analyte

9.1.2.1.2.2. QN_{istd} = response of the quantitative ion of the internal standard amount

9.2. Calibration

9.2.1. A linear regression resulting from the 6 standards is used to quantitate the analytes in the load. The area of the analyte divided by the area of the internal standard is used in the resulting formula of the calibration curve.

9.3. Dilution Factor

9.3.1. $D = \text{Total volume} / \text{Sample volume}$

9.4. Multiplier for homogenates, dilutions, and non-standard volumes

9.4.1. $M = (V_{\text{curve}} / V_{\text{samp}}) \times D$

9.4.1.1. M = Multiplier

9.4.1.2. D = dilution factor

9.4.1.3. V_{curve} = matrix volume of calibration curve

9.4.1.4. V_{samp} = matrix volume of specimen

9.5. Concentration

9.5.1. $C = (A / V) * M$

9.5.1.1. C = Concentration (ng/mL) of the analyte in the unknown case.

9.5.1.2. A = Amount of drug in sample

9.5.1.3. V = Volume of sample

9.5.1.4. M = Multiplier

9.6. Max/Min

9.6.1. Percent Difference = $((R_h / R_l) - 1) \times 100$

9.6.1.1. R_h = high result

9.6.1.2. R_l = low result

SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)

9.7. Average

9.7.1. Average = $(R_1 + R_2) / 2$

9.7.1.1. R_1 = first result

9.7.1.2. R_2 = second result

9.8. Qualifier Ion Ratios

9.8.1.1.1. Ratio 1 = QL_1/QN

9.8.2. QL_1 = response of the quantifying product ion

9.8.3. QN = response of the qualifying product ion

10. Quality Control

10.1. Acceptance criteria

10.1.1. Chromatogram

10.1.1.1. Peaks must be Gaussian shaped (symmetrical).

10.1.1.2. Peaks must not exhibit extreme fronting or tailing.

10.1.1.3. Peaks sharing parent/product ions must have baseline resolution.

10.1.1.4. The internal standard (ISTD) in each case should be inspected for evidence of signal enhancement and suppression. The area of the quantifying ion should not be less than 50% or more than 200% of the average ISTD of the calibrators.

10.1.1.5. Retention time must not deviate outside $\pm 3\%$ of target, based upon the retention time of the calibrators and controls.

10.1.2. Mass spectroscopy

10.1.2.1. The ion ratio of all samples must not be greater than $\pm 20\%$ of the target ratio, as determined by a mid-level calibrator (CAL 4).

10.1.2.2. Coelution of quantifying and qualifying ions must not be greater than 0.025 minutes.

10.1.3. Calibrators

10.1.3.1. Analytical curves must have a coefficient of determination (R^2) of 0.992 or greater.

SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)

10.1.3.2. Each calibrator, when calculated against the calibration curve, must not deviate outside $\pm 20\%$ of the target value.

10.1.3.3. Refer to “Calibration curve point exclusion guidelines” section of the QA/QC Manual.

10.1.4. Controls

10.1.4.1. Controls must calculate within $\pm 20\%$ of the target value.

10.1.5. Blanks

10.1.5.1. Blanks should not contain any target analyte signal with an internal standard response ratio greater than 10% that of the lowest calibrator for the same analyte.

10.1.6. Any deviation from the above criteria must be approved by a senior chemist.

11. Validation of Method

Parameter	Result
Bias	Blood - L: -4.11% H: 2.29% Liver - L: -1.22% H: -0.27%
Precision	Blood - L: 4.53% H: 4.55% Liver - L: 6.72% H: 3.11%
Calibration model	0.2 - 100 ng/mL - Linear (1/x)
Carryover	No carryover observed following injection of 2X highest calibrator (200ng/mL).
Interference Studies	No significant interfering signal from matrix, internal standard, common drugs of abuse (including metabolites), OTC drugs, and prescription medications was observed
Ionization/Suppression: (Not needed if IS coelutes within 0.05 min.)	Blood - Furanylfentanyl 125% fentanyl-d5 130% (5% diff). Liver - Furanylfentanyl 130% fentanyl-d5 126% (4% diff). All %RSD $\leq 20\%$ of target concentration.

SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)

LOD (Calculate: 3.3xSD Y-intercept/Mean of Slope)	0.04 ng/mL
LOQ (Set to lowest calibrator with acceptable Bias/Precision).	0.2 ng/mL
Dilution Integrity	NA - Specimens not routinely diluted
Processed Sample Stability - (re-analyze after 8 days)	Sample extracts stable for up to 8 days following extraction when refrigerated.
Recovery	Furanylfentanyl: Blood-86% Liver-66% Fentanyl-d5: Blood - 89% Liver - 72%

12. Reporting

12.1. The percent difference of duplicate analysis for an analyte must be less than or equal to 25% (see Max/Min in [Calculations](#) section).

12.2. Reporting of duplicate analysis should be done according to the table below:

Reporting Duplicates

- Dilution factors of 1 and 1

Scenario \ Dil	1	1	REPORT
A	In curve	In curve	Average
B	In curve	AQL or BQL	“In” value
C	In curve	ND *	Repeat
D	AQL/BQL	AQL/BQL	Less than/ Greater than
E	BQL	ND	ND

* ND = None Detected, due to IRC, S/N threshold, r.t., or other

12.2.1.1. In Curve = Measured concentration (pre-multiplier) falls within the calibration range

SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)

12.2.1.2. AQL = Measured concentration (pre-multiplier) falls Above Quantitation Limit

12.2.1.3. BQL = Measured concentration (pre-multiplier) falls Below Quantitation Limit

12.2.1.4. ND = None Detected

12.3. Averaging reportable values

12.3.1. Results for duplicate analysis (both falling within calibration curve) shall be truncated prior to averaging.

12.3.2. Enter calculated concentration for each specimen into toxlog.

12.4. Significant figures

12.4.1. Concentrations are truncated and reported with two significant figures in nanograms per milliliter (ng/mL).

13. Preparation of Load

13.1. The load paperwork and data is to be arranged in the following order:

13.1.1. Assignment sheet

13.1.2. Comments or note to file if applicable

13.1.3. Load summary

13.1.4. Specimen worklist

13.1.5. Chain of custody (Specimen)

13.1.6. Aliquot chain of custody

13.1.7. Standard and control worksheet

13.1.8. Sequence summaries/calibration reports – paper clipped

13.1.9. Calibrator data - paper clipped

13.1.10. Blank matrix data - paper clipped

13.1.11. Control data - paper clipped

13.1.12. Specimen data – stapled

SOP 416 – Furanylfentanyl Extraction using Protein Precipitation for Quantification by Liquid Chromatography/Electrospray Mass Spectrometry/ Mass Spectrometry (LC/MS/MS)

14. References

- 14.1. Helander, Anders, Matilda Bäckberg, and Olof Beck. "Intoxications Involving the Fentanyl Analogs Acetylfentanyl, 4-methoxybutyrfentanyl and Furanylfentanyl: Results from the Swedish STRIDA Project." *Clinical Toxicology* 54.4 (2016): 324-32.
- 14.2. Palamalai, Vikram, Kalen N. Olson, Julie Kloss, Owen Middleton, Kelly Mills, A. Quinn Strobl, Lindsey C. Thomas, and Fred S. Apple. "Superiority of Postmortem Liver Fentanyl Concentrations over Peripheral Blood Influenced by Postmortem Interval for Determination of Fentanyl Toxicity." *Clinical Biochemistry* 46.7-8 (2013): 598-602.
- 14.3. Olson, Kalen N., Kristin Luckenbill, Jonathan Thompson, Owen Middleton, Roberta Geiselhart, Kelly M. Mills, Julie Kloss, and Fred S. Apple. "Postmortem Redistribution of Fentanyl in Blood." *American Journal of Clinical Pathology Am J Clin Pathol* 133.3 (2010): 447-53.

Other suggested reading:

- 14.4. Chambers, Erin, Diane M. Wagrowski-Diehl, Ziling Lu, and Jeffrey R. Mazzeo. "Systematic and Comprehensive Strategy for Reducing Matrix Effects in LC/MS/MS Analyses." *Journal of Chromatography B* 852.1-2 (2007): 22-34.