

Deviation Request Form (DRF)

Directions: The Initiator will complete Sections A through C. Additional continuation pages can be included if necessary.

Initiator	Leslie Daugherty	Date	9/1/2020
A. Requested deviation applies to (Technical Procedure – include specific section):			
Technical Procedure for the use of the Cyanovac, Technical Procedure for Leuco Crystal Violet, Technical Procedure for Rhodamine 6G, Technical Procedure for the use of the Superglue Fuming Wands -3.0			
B. Requested deviation:			
Remove Ultra Lite ALS.			
C. Necessity for the deviation:			
No longer an available ALS.			
D. Technical review and Authorization (to be completed by the Quality Manager and/or Technical Leader)			
Comments(to include merits and impacts):			
Approved	<input checked="" type="checkbox"/>	Yes	<input type="checkbox"/> No
Duration		until processing procedure is approved	
Signature	Leslie Daugherty <small>Digitally signed by Leslie Daugherty DN: cn=Leslie Daugherty, e=ldaugherty@ncdoj.gov Reason: I am the author of this document Date: 2020-09-01 12:15:23 Foxit PhantomPDF Version: 9.3.0</small>		Date
		9/1/2020	
E. Quality Assurance Authorization (to be completed by the Quality Manager, Forensic Scientist Manager or designee)			
Acceptable within general QA guidelines and good laboratory practice?		<input checked="" type="checkbox"/>	Yes
		<input type="checkbox"/>	No
Significant negative impact to Crime Laboratory Quality System?		<input type="checkbox"/>	Yes
		<input checked="" type="checkbox"/>	No
Restrictions/limitations:			
<input checked="" type="checkbox"/>	Authorized	<input type="checkbox"/>	Rejected
Signature	Karen W. Morrow <small>Digitally signed by Karen W. Morrow Date: 2020.09.01 15:51:39 -04'00'</small>		Date
		9.1.2020	

Technical Procedure for the use of the Cyanovac

1.0 Purpose – This procedure is a non-porous development procedure.

2.0 Scope – This procedure is a step in the processing of non-porous evidence that may contain impressions that require developing/enhancing.

2.1 The Cyanovac polymerizes the latent impression using cyanoacrylate ester in a vacuum environment. The vacuum will eliminate background moisture and allow the Cyanoacrylate ester to attach to the components of the latent impression thereby eliminating the over-fuming that may occur with manual cyanoacrylate ester techniques. Numerous materials, including plastic bags, weapons, metals, and various other substrates, may be processed using the Cyanovac. Cyanoacrylate ester shall be used as a preliminary process when utilizing subsequent processing techniques. Fluorescent dye staining, in conjunction with laser examinations, is dependent on the proper use of cyanoacrylate ester fuming techniques.

3.0 Definitions

- **Alternate light source (ALS):** Any of the multiple forensic light sources readily available in the Latent Evidence Section including, but not limited to, the CrimeScope (CS), Mini CrimeScope (MCS), TracER Laser, and Ultra-Lite ALS.
- **Ambient light:** Light that is readily available in the office environment (i.e., natural light or light that emanates from an office lighting source).
- **Short tube:** Any of the Cyanovacs located in the Latent Evidence Section that are upright in nature and designated for use with shorter items of evidence.
- **Long tube:** Any of the Cyanovacs located in the Latent Evidence Section that are horizontal in nature and designed for use with longer items of evidence.
- **CE:** Cyanoacrylate ester, also known as super glue.
- **Cyanovac:** Any of the multiple Cyanovac units (long or short tube) that assist in the processing of non-porous items of evidence with cyanoacrylate ester in a vacuum environment.

4.0 Equipment, Materials and Reagents

4.1 Equipment and Materials

- Cyanovac (short or long tube)

4.2 Reagents

- Cyanoacrylate ester (Bottle/Vial)

5.0 Procedure

5.1 Items to be processed in the Cyanovac must first be examined visually in ambient lighting conditions and then with the use of an alternate light source (CS, MCS).

5.2 Remove the end cap from the Cyanovac chamber by releasing the elastic T-straps.

5.3 Insert item(s) of evidence into chamber by hand.

Note: Items may touch each other and the sides of the vessel without significant effect on processing. Large flexible items such as garbage bags must be unfolded to the degree possible, but need not be fully spread. Sealed items, such as zip top plastic bags or sealed plastic or metal containers must be opened in order to prevent rupture due to vacuum environment.

5.4 Place 5 to 10 drops of cyanoacrylate (super glue) in a foil dish and place into chamber.

5.5 Verify that the o-ring is present in the end cap. Ensure that it is not broken or dry-rotted.

5.6 Replace end cap on the chamber and fasten using elastic T-straps.

5.7 Verify that the chamber bleed valve is closed.

5.8 Turn the vacuum pump to ON using the switch.

5.9 Observe the vacuum gauge to ensure that the air in the chamber is being evacuated. Pump will automatically stop at the appropriate pressure set by the manufacturer.

Note: Pump will automatically restart if the vacuum pressure varies from optimum.

5.10 Allow the item(s) to remain under vacuum for 20 minutes.

Note: Some items may require a longer processing time; however, this period of time will not compromise the test value.

5.11 Turn vacuum pump to OFF using the switch.

5.12 Open the chamber bleed valve to equalize pressure.

5.13 Remove chamber end cap and remove cyanoacrylate source.

5.14 Leave item(s) in the chamber for approximately 10 minutes.

5.15 Remove item(s) and allow to rest for 24 hours.

Note: The additional 24 hour rest allows for full setting of cyanoacrylate.

5.16 Examine item(s) for developed latent prints using subsequent processing techniques.

5.17 Any developed latent prints must then be preserved using the method described in the Steps for Preserving Developed Impressions.

5.18 Standards and Controls – Forensic Scientists shall produce a self-made test print to be processed concurrently with items of evidence.

5.19 Calibration – See Cyanovac operating manual for further information on controls and specifications.

5.20 Sampling – N/A

5.21 Calculations – N/A

5.22 Uncertainty of Measurement - N/A

6.0 Limitations – Cyanovac is for use in the processing of non-porous evidence.

6.1 The cyanoacrylate fuming process is vital to subsequent treatment with fluorescent dyes and laser and/or alternate light source examinations (see Fluorescent Dyes and Laser/Alternate Light Sources).

7.0 **Safety** – Proper purging of the system is necessary as the fumes may cause irritation when in contact with the eyes or skin and may be harmful if inhaled or ingested. Protective goggles, gloves, and apron/lab coat shall be worn during processing. Additionally, cyanoacrylate ester is an adhesive/glue. Care shall be taken to avoid application to unintended surfaces.

8.0 References

Besonen, J.A. “Heat Acceleration of the Superglue Fuming Method for Development of Latent Fingerprints.” *Identification News*. (1983): 3 – 4.

Bessman, C.W., et al. “A Comparison of Cyanoacrylate Fuming in a Vacuum Cabinet to a Humidity Fuming Chamber.” *Journal of Forensic Identification*. Vol. 55, 1: 10 – 35 (2005).

Cummings, H., M. Hollars and T. Trozzi. “Getting the Most from Cyanoacrylate Dyes.” *Journal of Forensic Identification*. Vol. 43, 1: 37-43 (1993).

Day, K.J. and W. Bowker. “Enhancement of Cyanoacrylate Developed Latent Prints Using Nile Red.” *Journal of Forensic Identification*. Vol. 46, 2: 183-187 (1996).

Deobald, G.W. “The Effect of Cyanoacrylate Fuming on Firearms Examinations.” *Identification Canada*. (1992): 4-13.

Fallano, J.F. “Alternatives to Alternate Light Sources: How to Achieve a Greater Print Yield with Cyanoacrylate Fuming.” *Journal of Forensic Identification*. Vol. 42, 2: 91-95 (1992).

Fertgus, R.E. “Latent Print Destruction and Superglue Stabilization.” Florida Division of the International Association for Identification. (1993): 7.

Froude Jr., J.H. “The Super Glue Fuming Wand: A Preliminary Evaluation.” *Journal of Forensic Identification*. Vol. 46, 1: 19-31 (1996).

Howorka, H., and K. Kretschmer. “Experimental Study of Using Cyanoacrylate Ester Vapor for Developing Latent Fingerprints.” *Forensic Science International*. Vol. 46: 31-32 (1990).

Kendall, F.G., and B.W. Rehn. “Rapid Method of Super Glue Fuming for the Development of Latent Fingerprints.” *Identification News*. (June 1982): 3-4.

Kendall, F.G. “Superglue Fuming for the Development of Latent Fingerprints.” *Identification News*. (May 1982): 3-5.

Kent, T., ed. *Manual of Fingerprint Development Techniques: A Guide to the Selection and Use of Processing for the Development of Latent Fingerprints*. Police Scientific Development Branch, London (July 1992).

King, W.R. "The Effects of Differential Cyanoacrylate Fuming Times on the Development of Fingerprints on Skin." *Journal of Forensic Identification*. Vol. 59, 5: 537 – 544 (2009).

Kobus, H.J., R.N. Warren, and M. Stoilovic. "Two Simple Staining Procedures Which Improve the Contrast and Ridge Detail of Fingerprints Developed with "Super Glue" (Cyanoacrylate Ester)." *Forensic Science International*. Vol. 23: 233-240 (1983).

Lee, H.C. "Methods of Latent Print Development." *Proceedings of the International Forensic Symposium on Latent Prints*. (July 1987): 15 – 24.

Lennard, C.J. and P.A. Margot. "Sequencing of Reagents for the Improved Visualization of Latent Fingerprints." *Proceedings of the International Forensic Symposium on Latent Prints*. (July 1987): 141-142.

Llewellyn Jr., P.E. and L.S. Dinkins. "A New Use for an Old Friend." *Journal of Forensic Identification*. Vol. 45, 5: 498-503 (1995).

Manual of Fingerprint Development Techniques: A Guide to the Selection and Use of Processes for the Development of Latent Fingerprints. Scientific Research and Development Branch, London (1986).

Mazzella, W.D. and C.J. Lennard. "An Additional Study of Cyanoacrylate Stains." *Journal of Forensic Identification*. Vol. 45, 1: 5-18 (1995).

Menzel, E. R., et al. "Laser Detection of Latent Fingerprints: Treatment with Glue Containing Cyanoacrylate Ester." *Journal of Forensic Sciences*. Vol. 28, 2: 307-317 (April 1983).

Mock, J.P. "Cyanoacrylates and Heat – A Word of Caution." *The Identification Section*. Vol. 3, 3 (June 1985).

Olsen, R., and M.F. Shonberger. "A Slow-Reacting Catalyst for Cyanoacrylate Fuming." *Journal of Forensic Identification*. Vol. 45, 6: 651-653 (1995).

Sahs, P. T. and R.J. Wojcik. "Moisture Catalyst for Cyanoacrylate Fuming." *Identification News*. (September 1984): 9.

US Department of Justice. *Chemical Formulas and Processing Guide for Developing Latent Prints*. FBI Laboratory Division, Latent Fingerprint Section (1994).

Watkin, J.E., et al. "Cyanoacrylate Fuming of Latent Prints: Vacuum Versus Heat/Humidity." *Journal of Forensic Identification*. Vol. 44, 5: 545-556 (1994).

Weaver, D.E., and E. J. Clary. *A One Step Fluorescent Cyanoacrylate Fingerprint Development Technology*. State of Alaska Scientific Crime Detection Laboratory Research Team. 1-12.

Wilkinson, D.A., and A.H. Misner. "A Comparison of Thenoyl Europium Chelate with Ardrex and Rhodamine 6G for the Fluorescent Detection of Cyanoacrylate Prints." *Journal of Forensic Identification*. Vol. 44, 4: 387-406 (1994).

9.0 Records – N/A

10.0 Attachments – N/A

Revision History		
Effective Date	Version Number	Reason
09/17/2012	1	Original Document
10/31/2013	2	Added issuing authority to header
03/30/2017	3	Header Update – Removed Digital reference.
01/19/2018	4	Updated issuing authority in header 5.1 & 5.18 - Moved requirement for test print to “Standards and Controls.”
02/01/2019	5	Used the term “Cyanoacrylate Ester” instead of “Cyanoacrylate” throughout document 3.0: Added abbreviation for ALS, CrimeScope and Mini CrimeScope Changed number references to be numerals only instead of spelled out throughout document