Greiss Paper Solution

I. Required Items:

	<u>For 1000 mL</u>	<u>For 2500 mL</u>	<u>For 4000 mL</u>
Sulfanilic Acid	38.5 gr	96.25 gr	154.0 gr
Deionized water (diH ₂ O)	500 mL	1250 mL	2000 mL
Naphthol	21.5 gr	53.75 gr	86.0 gr
Methanol	500 mL	1250 mL	2000 mL

II. Instructions For Preparation:

Dissolve sulfanilic acid in diH_2O . Dissolve naphthol in methanol. Combine equal amounts of the sulfanilic acid and naphthol solutions.

III. Storage Conditions:

Store at room temperature.

IV. Expiration Date:

Expires one (1) year after preparation.

Vendor/Lot# Sulfanilic Acid	Vender/Lot# diH ₂ O	Vendor/Lot# Naphthol	Vendor/Lot# Methanol	Amount Made	Date Prepared	Preparer	Expiration Date

Version 2

Effective Date: 06/09/2014