Training Procedure for Preliminary Color Tests

Version 6

Effective Date: 08/17/2018

- **Purpose** A chemical color test is a preliminary screening test in which a sample of an unknown substance is added to a chemical reagent in order to produce a color change. The color change that occurs is indicative of a class or type of drug that may or may not be present. Knowledge of the structure of controlled substances can often be used to predict the results of color tests. In this section, the trainee will learn how color test reagents are used to classify and narrow the possible controlled substances present in forensic drug samples.
- **2.0 Scope** This procedure applies to trainees in the Drug Chemistry Sections of the State Crime Laboratory.
- 3.0 **Definitions** N/A
- 4.0 Procedure

4.1 Objectives

- **4.1.1** Review the Drug Chemistry Technical Procedure for Preliminary Color Tests and the Drug Chemistry Technical Procedure for Drug Analysis as it pertains to preliminary color tests.
- **4.1.2** Review the lab-wide Procedure for the Use of Forensic Advantage (FA) and complete an overview of FA) with an FA administrator within the drug chemistry section. This will include (but is not limited to) case submission details, data entry, and drug chemistry worksheet information.
- **4.1.3** Be able to document a newly prepared reagent and perform quality control procedures as listed in the Drug Chemistry Administrative Policy for Quality Assurance. Documentation will be performed in the resource manager portion of FA.
- **4.1.4** Be able to presumptively identify/classify what controlled/non-controlled substances may be present based on the results of performing the color tests.
- **4.1.5** Successfully complete a practice exercise (separately given in handout) and a written exam.

4.2 Study Ouestions

- **4.2.1** Define the term reagent.
- **4.2.2** What comprises a positive reaction for the Marquis reagent? What type of reaction is occurring here?
- **4.2.3** For the sulfuric acid based reagents, how does the forensic scientist know that the color change is not due to dyes? For liquid submissions, how does the forensic scientist know that the color change is not due to solvent interaction?
- **4.2.4** What functional group is generally associated with the Cobalt Thiocyanate reagent turning blue?

4.2.5 Why do cocaine base and cocaine hydrochloride react to the Cobalt Thiocyanate reagent with different intensities?

Version 6

Effective Date: 08/17/2018

- **4.2.6** What types of compounds give a positive reaction to pDMAB? Give some examples.
- **4.2.7** The Koppanyi and the Zwikker reagent generally test for what class of drugs?
- **4.2.8** The Duquenois-Levine reagent is used to test for what type of substance?
- **4.2.9** What causes the color change in the Duquenois-Levine reagent?
- **4.2.10** Are there any other substances (not controlled by CSA) that cause a positive reaction with the Duquenois-Levine reagent?
- **4.2.11** Explain why a color test cannot be used to identify a substance definitively.
- **4.2.12** How can you ensure that a color test reagent is actually working since some substances do not cause a color change? How does a Forensic Scientist avoid false positives and false negatives?
- **4.2.13** Why should unknowns be added in very small amounts (ex: what would happen if potassium cyanide was added to the spot well with Marquis Reagent?
- **4.2.14** How are primary and secondary reference materials obtained and documented for casework comparison purposes?

4.3 Practical/Laboratory Exercises

- **4.3.1** Prepare, label and properly document the following color test reagents: (See the Drug Chemistry Administrative Policy for Quality Assurance for assistance.)
 - Marquis (Mecke and Froehde optional)
 - Cobalt thiocyanate
 - Secondary amine
 - Koppanyi
 - Zwikker
 - p-Dimethylaminobenzaldehyde (PDMAB)
 - Ferric Chloride
 - Dugenois-Levine
 - Liebermann's
 - Chen-Kao
- **4.3.2** Test a set of known training standards provided to you by the Drug Chemistry Training Coordinator or designee and record the results. (Availability of standards may necessitate substitution of alternate substances. This is just a general guideline of common substances and subsequent presumptive testing.)

(**Note:** Standards listed for testing with the Marquis reagent should be also be tested with concentrated sulfuric acid so the Forensic Scientist Trainee can see how a strong acid may react with these compounds. Also, if not noted specifically, please test compounds

Version 6 Effective Date: 08/17/2018

with both secondary amine #1 and #2)

4.3.2.1 Standard Preliminary Color Test

Amphetamine HCl Marquis, cobalt, secondary amine, Chen-Kao Methamphetamine HCl Marquis, cobalt, secondary amine, Chen-Kao

Dimethylsulfone (DMSO₂₎ Marquis, cobalt

Phentermine Marquis, cobalt, secondary amine, Chen-Kao

Pseudoephedrine HCl Marquis, cobalt, Chen-Kao

Aspirin Marquis, cobalt

Cocaine HCl Marquis, cobalt, PDMAB, secondary amine #1
Cocaine base Marquis, cobalt, PDMAB, secondary amine #1

Diphenhydramine Marquis, cobalt

Fentanyl (or derivative of) Marquis, secondary amine #1

Guaifenesin Marquis, cobalt

Heroin base Marquis, cobalt, secondary amine #1
Morphine base Marquis, cobalt, secondary amine #1
Methadone HCl Marquis, cobalt, secondary amine #1
MDA Marquis, cobalt, secondary amine, Chen-Kao
MDMA Marquis, cobalt, secondary amine, Chen-Kao

Inositol Marquis, cobalt
Mannitol Marquis, cobalt
Sodium bicarbonate Marquis, cobalt
Toluene Marquis, cobalt
GBL (if available) Ferric Chloride, cobalt
GHB Ferric Chloride, cobalt

BZP Marquis, cobalt, secondary amine #1 TFMPP Marquis, cobalt, secondary amine #1

Benzocaine Marquis, cobalt, PDMAB, secondary amine #1

Lidocaine HCl Marquis, cobalt, PDMAB Procaine HCl Marquis, cobalt, PDMAB

LSD PDMAB, check for fluorescence with long wave UV light

LAMPA PDMAB

Amobarbital Marquis, Zwikker, Koppanyi Marquis, Zwikker, Koppanyi Pentobarbital Secobarbital Marquis, Zwikker, Koppanyi Hashish Duquenois-Levine, Liebermanns Duquenois-Levine, Liebermanns Marijuana Patchouli oil Duquenois-Levine, Liebermanns Tobacco Duquenois-Levine, Liebermanns Synthetic cannabinoid (if available) Duquenois-Levine, Liebermanns 3,4-Methylenedioxypyrovalerone Marquis, Liebermanns, Chen-Kao

(MDPV)

Cathinone (or substituted) Marquis, secondary amine, Liebermanns, Chen-Kao

4.3.3 Be able to give an oral explanation of color testing in a manner that would be applicable to a jury trial.

5.0 Required Reading

O'Neal, C. L. et al. "Validation of twelve chemical spot tests for the detection of drugs of abuse." *Forensic Science International*, 109 (2000): 189-201.

Toole, K.E. et al. "Color Tests for the Preliminary Identification of Methcathinone and Analogues of Methcathinone." *Microgram Journal*, Volume 9, Number 1 (Jan. 1999): 27-32.

Tsumura, Y. et al. "False positives and false negatives with a cocaine-specific field test and modification

Version 6

Effective Date: 08/17/2018

6.0 References

Clarke, E.G.C. Isolation and Identification of Drugs. London: The Pharmaceutical Press, 1969.

of test protocol to reduce false decision." Forensic Science International, 155 (2005): 158-164.

M. Windholz, ed. The Merck Index. 9th Ed. Rahway: Merck & Co., 1976.

O'Neal, C. L. et al. "Validation of twelve chemical spot tests for the detection of drugs of abuse." *Forensic Science International*, 109 (2000): 189-201.

Pitt, Coltin G. *Investigation of problems Involving Chemical Analysis of The Duquenois Levine Color Reaction*. Research Triangle Institute: RTI Project 31U-628. November 1971.

Poyner, B. et al. "Presumptive Color Test for Synthetic Cannabinoids Containing an Indole Substructure." *Journal of the Clandestine Laboratory Investigating Chemists Association*. Volume 22 Number 4 (Oct. 2012): 27-31.

Toole, K.E. et al. "Color Tests for the Preliminary Identification of Methcathinone and Analogues of Methcathinone." *Microgram Journal*, Volume 9, Number 1 (Jan. 1999): 27-32.

Tsumura, Y. et al. "False positives and false negatives with a cocaine-specific field test and modification of test protocol to reduce false decision." *Forensic Science International*, 155 (2005): 158-164.

7.0 Records

- Drug Chemistry Training Checklist
- Section Completion Summary

8.0 Attachments - N/A

Revision History			
Effective Date	Version Number	Reason	
09/17/2012	1	Original Document	
05/03/2013	2	Scope – Reworded to match other ISO documents.	
		4.1.1 – Preliminary added to technical procedure name; added Drug Chemistry Technical Procedure for Drug Analysis, as it pertains to preliminary color tests.	
		4.1.2 – Name of Administrative Policy for Quality Assurance updated.	
		4.3.1 – Name of Administrative Policy for Quality Assurance updated. Bullet list – (Mecke and Froehde optional) added with Marquis. Added Ferric Chloride and Liebermann's to list.	
		4.3.2 – Added (Availability of standards may necessitate substitution of alternate substances.) Added note reference use of	

		sulfuric acid for standards tested with Marquis reagent. 4.3.2.1 bullet list – Added GHB and GBL with Ferric Chloride and cobalt tests, added Liebermann's to Synthetic cannabinoid, and added MDPV with Liebermann's. Added check for fluorescence with long wave UV light to LSD. References – Added Poyner and Toole references for Liebermann's Reagent.
12/06/2013	3	Added issuing authority to header
10/19/2015	4	Header – Revised issuing authority 4.1.3 – Added new objective reference primary and secondary standards
		4.2.3 – Added new study question reference sulfuric acid based reagents and dyes
		4.2.11, 4.2.13 – Edited study questions
		4.2.19 – Added study question reference primary and secondary reference standards.
		Required Reading – Added new section
		References – Added additional references
08/30/2016	5	4.2.7 – updated coloring to yellow 4.2.10 – fix typo
08/17/2018	6	4.1.2 – added objective; objective numbers changed 4.1.3, 4.1.5, 4.1.6 – updated wording 4.1.4 - removed 4.2 – updated study questions 4.2.12, 4.2.14, 4.2.15 – removed these questions 4.3.1 – added Chen-Kao 4.3.2 – updated with new testing/wording 4.3.3 - added 7.0 – removed reagent log (paper copy)
		_
		_

Version 6

Effective Date: 08/17/2018