CHARLOTTE-MECKLENBURG POLICE DEPARTMENT CRIME LAB CHEMISTRY SECTION

STANDARD OPERATING PROCEDURES

Controlled Substance Analysis Standard Operating Procedures

Section 1Procedures:

- 1. General Evidence Procedures
 - 1.1 Introduction
 - 1.2 Evidence Acceptance
 - 1.3 Evidence Handling
- 2. Analysis
 - 2.1 Introduction
 - 2.2 Analytical Standards
 - 2.3 Analytical Schemes
 - 2.4 Preliminary Testing
 - 2.5 Confirmatory Testing
- 3. Measurement Practices
 - 3.1 Sampling
 - 3.2Weighing
 - 3.3 Uncertainty of Measurement
- 4. Reporting and Documentation Procedures
 - 4.1 Case Files
 - 4.2 Case Notes, Worksheets, and Photographs
 - 4.3 Reports
- 5. Quality Assurance
 - 5.1 General
 - 5.2 Reagents
 - 5.3 Standards
 - 5.4 Blanks
 - 5.5 Color Tests
 - 5.6 Balances
 - 5.7 Validation of Methods/Procedures
 - 5.8 GS/MS
 - 5.9 FT-IR
 - 5.10 Gases
 - 5.11 Peer Review

Section 2 Appendices:

- A. Abbreviations Commonly Used
- B. Color Tests
- C. Unit Conversions
- D. Hypergeometric Sampling Table
- E. NC Trafficking Limits
- F. Commonly Used Instrument Methods
- G. Special Procedures
- H. Equipment
- I. References

Fire Debris Analysis Standard Operating Procedures

Section 1 Procedures

- 1 General Evidence Procedures
 - 1.1 Introduction
 - 1.2 Evidence Acceptance
 - 1.3 Evidence Handling
- 2 Scheme of Analysis
- 2.1 Introduction
- 2.2 General Scheme of Analysis
- 3 Sample Preparation/Techniques of Analysis
 - 3.1 Introduction
 - 3.2 Sample Preparation
- 4 Techniques of Analysis
 - 4.1 Introduction
 - 4.2 GC/MS
 - 4.3 Flammability Determination
 - 4.4 Miscibility Test
- 5 Data Interpretation and Classification of Ignitable Liquids
 - 5.1 Introduction
 - 5.2 Data Interpretation
 - 5.3 Extracted Ion Chromatograms
 - 5.4 Classification of Ignitable Liquids
- 6 Reporting and Documentation Procedures
 - 6.1 Introduction
 - 6.2 Case Files
 - 6.3 Case Notes, Worksheets, and Photographs
 - 6.4 Reports
- 7 Quality Assurance
 - 7.1 General
 - 7.2 Reagents
- 7.3 Standards
- 7.4 Blanks
- 7.5 System positive/Negative controls
- 7.6 Validation of Methods/Procedures
- 7.7 GC/MS
- 7.8 Gases
- 7.9 Supplies
- 7.10 Peer Review

Section 2 Appendices

- Appendix A: Abbreviations Commonly Used
- Appendix B: Commonly used Instrument Methods
- Appendix C: References
- Appendix D: Fire Debris Macro

Blood Analysis for Ethyl Alcohol by Headspace Gas Chromatography

- 1. Principle
- 2. Analytes Detected
- 3. Materials
 - 3.1. Control and Calibrators
- 4. Procedure
- 5. Apparatus and Materials
- 6. Headspace/Agilent 7890 GC/FID
- 7. Calculations
- 8. Quality Control
- 9. Linearity
- 10. Reporting
- 11. Uncertainty of Measurement
 - 11.1. General
 - 11.2. Measurement of Uncertainty
- 12. References

Chemistry Training Program
Fire Debris Training Program
Blood Alcohol Analysis Training Manual

Official Abbreviations

Document History

Version #	Issue Date	History
1.	03/01/00	Original Issue
2.	08/14/00	Change of GC calibration and addition of Ethylene Gycol method
3.	02/05/01	Validation of 5973 of MSD, UV, and balance cal., re-org. of
		document, expansion safety forms, addition of GCMS library SOP, val of
		5970 MSD cal., val of confirmatory GC idents., val of methods for 5973
		MSD
4.	02/23/01	changes in GHB method and LSD method
5.	09/05/01	changes in UV calibration procedure
6.	11/26/01	changes in FTIR and MS calibration procedures
7.	02/21/02	changes in MS calibration, Thermometer cal., Records keeping, UV
		classification (ancillary), BA evidence storage, typo corrections, GC
		calibration forms .
8.	11-26-02	BA reporting, safety checks assignments, weighing issues, rep sample
		issues.
9.	12-03-03	BA process, Mass Spec calib, GC calib., evidence maint log elimination,
		BA prof test, reprt/lab file
10.	11-23-04	change MS calibration procedure.
11.	02-25-05	method for sulfates, carbonates, posphates.
12.	08-19-05	new FTIR calibrations, Rt window changes, BA std change, admin
		changes
13.	04-12-06	delete in house arson prof test
14.	09/07/06	fixed table that was cut off in last update - no method changes
15.	05/31/07	changes to Analysis functions (blood alcohol) and operational functions
		(instrument analysis and testing format and controlled substance
		analysis), updated abbreviations page
16.	04/10/08	Changed technical review requirements for BA and Fire Debris Analysis
17.	12/09/08	* Major Revision of the entire SOP
18.	05/07/09	Delete Monthly Balance Calibrations and Thermometer Calibrations
19.	02/08/10	*Major Revision of the entire SOP and the addition of Fire Debris
		Analysis
20.	04/07/10	Updated abbreviations page, deleted outdated sections of GC-FID,
		Changes to peer review in quality section.
21.	07/20/11	*Major Revision of the entire SOP
22.	09/08/11	Fire Debris SOP:changes to Evidence Acceptance and Quality
		Assurance, Abbreviations List removed and made separate document
		Controlled Substance SOP: changes to Evidence Acceptance, Sampling,
		Reports, Quality Assurance General, QA Color Tests, QA Balances, QA
		GC/MS, QA FT-IR, Abbreviations List removed and made separate
		Document
		Alcohol Analysis SOP: changes Evidence Receiving and Handling,
		Sample Treatment, Evaluation of Quality Control Data, Documentation,
		Recording and Reporting
23.	12/15/11	No substantive changes made. However, References were added and
		the spacing, margins and font styles were fixed and made more uniform.
24	2/6/12	5.13 removed from the ISO Drug SOP.
25	2/27/12	Alcohol Analysis by Direct Inject had "should" changed to "may" for
		Calibration #2 and Fire Debris SOP had 6.9.3 changed to 10% technical
		review.
26	2/29/12	Blood Analysis for Ethyl Alcohol by Headspace Gas Chromatography
		SOP added and a Fire Debris Training Program was added.

27	3/22/12	Added headspace Load list form, added Multivolatile standard and
20	2/20/42	control. Changed "lid" to "cap"; added "Standby" on sequence run table.
28	3/30/12	Added the Blood Alcohol Training Program.
29	6/04/12	Drug Chemistry - Added Liebermans test (color tests), changed Marquis expiration to six months, removed 3.1.11; BAC- added sample treatment, minimum standards used for calibration, chromatography acceptance criteria - to Direct Inject and Headspace methods; and Fire debris had major or complete changes to 1, 3, 4, 5, 6 and 7.
30	8/29/12	Fire Debris – changes in 4.2.6, 5.3 (added a table), 7.7.2, Appendix B; Appendix D –added Fire Debris Macros
31	12/11/12	Measurement practices (3.1.3 – 3.2.17), Secondary Standards (5.3.7); MS tunes (5.8.4.8); added Simon's test – color reagents (Appendix B); Added steroid method- appendix F; added new inst. Galileo GC/MS (appendix H).
32	02/01/13	Took off 5.8.4.7; FTIR val-Q day of use or after maintenance (5.9.3); run caffeine weekly and mixed standard monthly (5.8.4).
33	12/09/13	Drug analysis UoM calculated, Balance performance/validation changed 5.6
34	01/08/14	Added transportation and handling of weights. (5.6.1)
35	04/07/14	Drug ID (Fire Debris)- Added PLIMS electronic Request 1.1.0(1.2.1); Evidence transfer using PLIMS 1.2.1(1.3); Lab number 1.2.4(1.3); Electronic case Reports/records 4.1 (6.2); Peer Review 5.11 (7.10) Fire Debris – CMPDF-(D1.MAC, D2.MAC, and D3.MAC) changed SAMPLE to Sample Name BAC – PLIMS electronic request/transfer; instrument batch creation and Data processing.
36	01/06/14	Removed Direct Inject from TOC; Changed FTIR Val-Q to Val-Pro per new software 5.9; added PLIMS to BAC, Fire debris and Drug analyst training manuals.
37	04/14/15	BAC headspace SOP: Added a new control (0.08mg/dl), moved to a 4-point calibration, serum/plasma – whole blood calculation.

Approval

Director	Matthew C. Mathis	Date:	
<u>Issuance</u>			
Section Administrator	Lillian Ngong	Date:	
Criminalist	Ann M. Charlesworth	Date:	
Criminalist	Andrew P. Oprysko	Date:	