

25% Nitric Acid Solution

I. Required Items:

	<u>For 300 mL</u>
Nitric Acid	75 mL
Deionized water (diH ₂ O)	225 mL

II. Instructions For Preparation:

Add nitric acid to diH₂O and mix.

III. Storage Conditions:

Store at room temperature in the top tray of the serial number kit or appropriate storage cabinet.

IV. Expiration Date:

Expires one (1) year after preparation.

Vendor/Lot# Nitric Acid	Vender/Lot# diH ₂ O	Amount Made	Date Prepared	Preparer	Expiration Date