

2N Sodium Hydroxide Solution, Method 1

I. Required Items:

	<u>For 100 mL</u>	<u>For 500 mL</u>
10N Sodium Hydroxide	20 mL	100 mL
Deionized water (diH ₂ O)	80 mL	400 mL

II. Instructions For Preparation:

Combine 10N sodium hydroxide with diH₂O. The volume prepared should be only the necessary volume for appropriate waste disposal.

III. Storage Conditions:

N/A – Only enough solution is made for immediate use.

IV. Expiration Date:

N/A – Only enough solution is made for immediate use.

Vendor/Lot# 10N Sodium Hydroxide	Vender/Lot# diH ₂ O	Amount Made	Date Prepared	Preparer	Expiration Date