

2N Sodium Hydroxide Solution, Method 2

I. Required Items:

Sodium Hydroxide
Deionized water (diH₂O)

For 100 mL
123.46 gr
100 mL

II. Instructions For Preparation:

Dissolve sodium hydroxide in diH₂O until the pellets are no longer visible. The volume prepared should be only the necessary volume for appropriate waste disposal.

III. Storage Conditions:

N/A – Only enough solution is made for immediate use.

IV. Expiration Date:

N/A – Only enough solution is made for immediate use.

Vendor/Lot# Sodium Hydroxide	Vender/Lot# diH ₂ O	Amount Made	Date Prepared	Preparer	Expiration Date