

5% Cartridge Case Cleaning Solution

I. Required Items:

	<u>For 100 mL</u>	<u>For 500 mL</u>
Birchwood Casey	5 mL	25 mL
Deionized water (diH ₂ O)	95 mL	475 mL

II. Instructions For Preparation:

Combine Birchwood Casey Cartridge Case Cleaner with diH₂O.

III. Storage Conditions:

Store at room temperature.

IV. Expiration Date:

Expires one (1) year after preparation.

Vendor/Lot# Cartridge Case Cleaner	Vender/Lot# diH ₂ O	Amount Made	Date Prepared	Preparer	Expiration Date