



# Specialty Tubes CONT'D

BD Vacutainer® Specialty Tubes

| Reference #   | Description  | Tube Size (mm)   | Draw Volume (mL) | Closure |        | Label type |
|---|--|------------------|------------------|---------|--------|------------|
|   |  |                  |                  | Type    | Colour |            |
| <b>Specialty - Sedimentation Rate Determination (SDR)</b> |  |                  |                  |         |        |            |
| 366676  | G Tube with Buffered Citrate .105M, CE Marked  | 8x100            | 1.8              |         |        | P          |
| 369741  | G Tube with Buffered Citrate .105M   | 13x75            | 2.4              |         |        | P          |
| 366465  | G Tube with Buffered Citrate .129M   | 13x100           | 4.0              |         |        | P          |
| 366065*   | G Tube with Buffered Citrate .105M, CE Marked  | 10x120           | 5.0              |         |        | P          |
| 366666*   | G Tube with Buffered Citrate .105M, CE Marked  | 10x120           | 5.0              |         |        | P          |
| <b>Specialty - Sodium Polyanethanesulfonate (SPS)</b>     |  |                  |                  |         |        |            |
| 366404  | G Whole Blood Tube with Additive 0.68ml  | 10.25x82         | 3.32             |         |        | P          |
| 364960  | G Whole Blood Tube with Additive 1.7ml   | 16x100           | 8.3              |         |        | ST         |
| <b>Specialty - Thrombin</b>                               |  |                  |                  |         |        |            |
| 367755  | G Serum tube, silicone-coated interior with Thrombin 1.4 NIH units                           | 13x100           | 7.0              |         |        | P          |
| 366525  | G Serum tube, silicone-coated interior with Thrombin 2.0 NIH units                           | 16x100           | 10.0             |         |        | P          |
| <b>Specialty - Miscellaneous</b>                          |  |                  |                  |         |        |            |
| 367855  | Tube with K2EDTA (spray-dried) for Lead Testing  | 13x75            | 3.0              |         |        | P          |
| 366401  | G Two tubes packaged together in a sterile pouch - Serum tube, silicone coated - K3EDTA tube | 16x100<br>13x100 | 10.0<br>7.0      | <br>    | <br>   | P<br>P     |
| 367001  | G Tube with 100mg Sod Fluoride/20mg Pot Oxalate, for Blood Alcohol Testing                   | 16x100           | 10.0             |         |        | P          |

\* For use with BD Seditrainer™ Rack (366016)

Specialty Tubes Continued on Next Page



3. Expiration Date

2. Lot Number

1. Reference/Catalogue Number: 367001

The vacutainer pictured here has the same ref number as the items that were recalled, but has a different lot number and expiration date from the affected lot. A vacutainer from the affected lot will have: 1) Reference/catalogue number: 367001, 2) Lot number: 8187663, and 3) Expiration date: 7/31/2020.