CHARLOTTE - MECKLENBURG POLICE DEPARTMENT

Crime Laboratory – Biology Section 1 M Dithiothreitol (DTT)



1M Dithiothreitol (DTT)

 Required items 	Э.
------------------------------------	----

DTT 0.77g DI water 5 ml

2. Instructions for Preparation:

Mix the DTT and DI water.

3. Storage Conditions:

Aliquot and freeze at -20°C. Do not freeze again once defrosted.

4. Expiration Date:

When the individual component expires, or twelve months after preparation, whichever comes first.

5. QC needed.

Lot # DTT Exp. Date	Amount Made	Date Made	Analyst	Sol'n Exp. Date	QC Date