



Genetic Analyzer AB 3130 Spectral

3130 # _____

For instructions, refer to –

- Getting Started Guide: Chapter 3 – Performing a Spectral Calibration
- Product Insert for Matrix Kit DS-33 Dye Set G5

Date Performed: _____

Name of Analyst: _____

Pass Fail (Check One)

Spectral Calibration: Required once a year unless needed before due to the following:

Check Reason:

- Annual
- A new dye set is used on the instrument
- The capillary length or polymer type is changed
- After the laser or CCD camera has been realigned/replaced by a service engineer
- When see a decrease in spectral separation (pull-up and or/pull down peaks) in the raw or analyzed data

Chief Criminalist/DNA Team Leader (or designee) Approval: _____

Comments: