



**Kernechtrot Stain (NFR)**

**1. Required items:**

Aluminum Sulfate	5.0 g
Nuclear Fast Red	0.1 g
DI water	100 ml

**2. Instructions for Preparation:**

Dissolve the Aluminum Sulfate in the hot DI water and add the Nuclear Fast Red. Stir and allow to cool, then filter.

**3. Storage Conditions:**

Room Temperature

**4. Expiration Date:**

When the individual components expire, or twelve months after preparation, whichever comes first.

Lot # Aluminum Sulfate/ Exp. Date	Lot # Nuclear Fast Red/Exp. Date	Amount Made	Date Made	Analyst	Solution Exp. Date