

**Proteinase K (13.5 mg/ml)****1. Required items:**

Pro K 100 mg
DI water 7.4 ml

2. Instructions for Preparation:

Add Pro K to the water. Aliquot small portions into tubes.

3. Storage Conditions:

Freeze at -20°C. Once thawed, do not re-freeze.

4. Expiration Date:

When the individual component expires, or twelve months after preparation, whichever comes first.

5. QC needed.

Lot # Pro K (Batch#) Exp. Date	Amt. Made	Date Made	Analyst	Sol'n Exp. Date	QC Date