

15% Acetic Acid Solution

I. Required Items:

	<u>For 1000 mL</u>	<u>For 4000 mL</u>
Acetic Acid (99.7+%, A.C.S. Reagent)	150 mL	600 mL
Deionized water (diH ₂ O)	850 mL	3400 mL

II. Instructions For Preparation:

Combine acetic acid with diH₂O.

III. Storage Conditions:

Store at room temperature.

IV. Expiration Date:

Expires six (6) months after preparation.

Vendor/Lot# Acetic Acid	Vender/Lot# diH ₂ O	Amount Made	Date Prepared	Preparer	Expiration Date