

5% Hydrochloric Acid Solution

I. Required Items:

	<u>For 1500 mL</u>	<u>For 2000 mL</u>	<u>For 4000 mL</u>
Hydrochloric Acid	75 mL	100 mL	200 mL
Deionized water (diH ₂ O)	1425 mL	1900 mL	3800 mL

II. Instructions For Preparation:

Combine hydrochloric acid with diH₂O.

III. Storage Conditions:

Store at room temperature.

IV. Expiration Date:

Expires six (6) months after preparation.

Vendor/Lot# Hydrochloric Acid	Vender/Lot# diH ₂ O	Amount Made	Date Prepared	Preparer	Expiration Date