

## **Ferric Chloride Solution**

### **I. Required Items:**

	<u>For 100 mL</u>	<u>For 300 mL</u>
Ferric Chloride	385.81 gr	1157.43 gr
Deionized water (diH <sub>2</sub> O)	100 mL	300 mL

### **II. Instructions For Preparation:**

Combine ferric chloride with diH<sub>2</sub>O and stir until the ferric chloride is no longer visible in the solution. The solution will appear as an orange color.

### **III. Storage Conditions:**

Store at room temperature in the top tray of the serial number kit or appropriate storage cabinet.

### **IV. Expiration Date:**

Expires one (1) year after preparation.

Vendor/Lot# Ferric Chloride	Vender/Lot# diH <sub>2</sub> O	Amount Made	Date Prepared	Preparer	Expiration Date