



Certificate of Quality Control

DNA IQ Kit

DNA IQ Reference Sample Kit

DTT

PFE Kit

PFE BTA Kit

ProK¹

QIAamp DNA Investigator Kit¹

¹(archived QIAamp blood extraction procedure will be used for QC)

Specifications:

QAM 8.2.4.5: A known human sample, preferably a NIST or NIST-traceable standard, will be extracted and quantified. To pass QC, at least 0.075ng/μl of DNA (0.015ng/μl for DNA IQ Reference Sample Kit) must be observed in the sample and no DNA observed in the reagent blank.

Lot # _____

Expiration Date _____

This is the expiration after testing, not the manufacturer expiration date.

For DNA IQ & QIAamp only:	For DNA IQ only:
Ethanol Lot # _____ Exp. _____	Propanol (2-) Lot # _____ Exp. _____

QC analyst _____

Date QC completed _____

Required records:

- Certificate of Quality Control
- DNA extraction worksheet
- QPCR results

Results of Testing: Pass Re-test (Choose One)

Results of Re-testing: Pass Fail Not Needed (Choose One)

Comments:

Approval of COQC Record
by Chief Criminalist, DNA
Team Leader, or designee: _____

Date* _____

*If passed, this is the date released