Ş

Picroindigocarmine Stain (PICS)

1. Required items:

Picric Acid	4.0 g
Indigocarmine	0.1 g
DI water	300 ml

2. Instructions for Preparation:

To DI water add Picric Acid and stir overnight. Add Indigocarmine and stir for at least one hour or until dissolved.

3. Storage Conditions:

Room Temperature

4. Expiration Date:

When the individual components expire, or twelve months after preparation, whichever comes first.

Lot # Indigocarmine/ Exp. Date	Amount Made	Date Made	Analyst	Solution Exp. Date
	Lot # Indigocarmine/ Exp. Date	Lot # Amount Made Indigocarmine/ Made Exp. Date Indicator Indigocarmine/ Indicator Indicator Indicator	Lot # Indigocarmine/ Exp. Date Amount Made Date Made Image: Date Image: Date Image: Date Image: Date Image: Date Image: Date <	Lot # Indigocarmine/ Exp. Date Amount Made Date Made Analyst Image: Constraint of the second of the